

ICS-2000 Ion Chromatography System



The ICS-2000 system is the first totally integrated and preconfigured Reagent-Free™ Ion Chromatograph designed to perform all types of electrolytically generated isocratic and gradient IC separations using conductivity detection. The ICS-2000 is available with a dual-piston pump, LCD touch-pad front panel, Reagent-Free eluent generation, thermally controlled conductivity cell, column heater, and optional vacuum degas. When coupled with AutoSuppression®, the ICS-2000 system provides high performance with unequalled ease of use. Automation provides full control and digital data collection from a PC using USB, high-speed communication protocol.

Versatile

- Performs all types of IC separations using conductivity detection.
- Integrated, preconfigured, factory plumbed, and tested for immediate productivity.
- Reagent-Free technology converts deionized water into high-purity eluents on-line.
- Streamlined design with small footprint occupies minimal bench space.
- An LCD touch-pad front panel provides clear identification of key operating parameters permitting at-instrument control and monitoring.
- Dual-piston pump design reduces pulsations, allowing high-sensitivity detection and excellent flow-rate accuracy and precision.

Reagent-Free IC

- Isocratic and gradient electrolytic eluent generation is capable of delivering hydroxide eluents for anion separations and methane-sulfonic acid eluents for cation separations.
- Automated eluent generation minimizes time, labor, operation costs, and eluent preparation errors.

Simple and Precise Control

- Built-in control for the SRS® and Atlas® electrolytic suppressors. AutoSuppression with electrolytic suppression eliminates the need to hand-prepare acid or base regenerants. Electrolytic suppression reduces background conductivity and provides high signal-to-noise ratios.



- Full control and digital data collection available with Windows® based Chromeleon® Chromatography Workstation Software using a USB high-speed communication protocol.
- Application templates preload all instrument parameters for fast and easy operation.
- Through Chromeleon software control—an electronic logbook provides monitoring of user-selectable operational parameters by creating virtual channels.

High Performance

- For improved reproducibility, the heated and thermostated high-performance conductivity detection cell permits measurements that are unaffected by temperature variation.


- Advanced single-range digital output with operating range to 3200 μS full scale. Alternate mode permits single-range analog signal output.
- Column heater provides day-to-day consistency, ensuring reproducibility and stability. Preheating of the eluent prior to the column ensures the column temperature set by the user. A transparent cover allows viewing of the column without temperature disruption.
- Optional built-in vacuum degas provides in-line degassing of eluents ensuring, reproducibility and protection of eluents from contamination and decomposition. Control of the degas operation can be automated to sense when degassing is required.
- Inert, nonmetallic PEEK components throughout the system ensure compatibility and metal contamination-free chromatography.

Convenient

- Versatile eluent organizer tray accommodates 1-, 2-, or 4-liter eluent bottles.
- Electrically actuated six-port Rheodyne PEEK injection valve.
- Ergonomically placed injection port for easy manual sampling.
- Eluent valve provides positive shut-off of eluent flow prior to the pump for easy servicing.
- Easy-access door to chromatography components.
- Leak detection and management allow fast response to system leaks.
- TTL controls for external pump, injection valve, range selection, and signal offset for stand-alone operation.

Key Features

- Reagent-Free IC
- LCD front panel control
- Dual-piston pump
- Column heater
- Electrolytic suppression
- Digital conductivity detection
- Vacuum degas (Option)
- USB connectivity, plug-n-play
- Optical leak detector
- Electronic logbook and trending through virtual channels



Reagent-Free IC produces consistent lab-to-lab eluent concentrations for highly reproducible retention times and peak areas. Results are the same day-to-day, system-to-system, and lab-to-lab.

All components are easily accessed through the front chromatography panel.

ICS-2000 IC SYSTEM SPECIFICATIONS

Analytical Pump and Hydraulics

Type:

Serial dual reciprocating pistons, microprocessor-controlled constant stroke, variable speed

Construction:

Chemically inert, metal-free PEEK pump heads and flow paths compatible with aqueous eluents of pH 0–14 and reversed-phase solvents

Control Mode:

Full control through front panel or through Chromeleon software; alternative control through TTL or relay closures

Maximum Operating Pressure:

35 MPa (5000 psi)

Pressure Ripple:

<1.0 % from 0.4 to 2.0 mL/min at 1000–3000 psi

Flow-Rate Range:

0.05–5.0 mL/min in 0.01 increments. Typical operating range is 0.4–2.0 mL/min

Flow Precision:

<0.2%

Flow Accuracy:

<1.0 % of set value or ± 2 μ L/min, whichever is greater

Piston Seal Wash:

Dual-pump head, wash can be continuous when connected to rinse solution supply

Pressure Alarm Limits:

Upper limit –35 MPa or 0–5000 psi in one unit (MPa or psi) increments; lower limit can be set up to one unit lower than upper limit

Vacuum Degas (Option):

Single channel, automatic

Conductivity Detector Electronics and Flow Cell

Type:

Microprocessor-controlled digital signal processor

Cell Drive:

8 kHz square wave

Linearity:

1% at 1 mS

Resolution:

0.1 nS

Full-Scale Output Ranges:

Digital signal range 0–3200 μ S

Analog signal range 0–3000 μ S

Electronic noise:

± 0.1 nS when background conductivity is 0–150 μ S

± 2 nS when background conductivity is 151–3200 μ S

Filter:

Rise times from 0 to 10 s, user selectable

Temperature Compensation:

Fixed at 1.7% per 1 °C at cell temperature

Temperature Range:

Ambient +7 °C to 55 °C

Temperature Stability:

≤ 0.01 °C

Cell Electrodes:

Passivated 316 stainless steel

Cell Body:

Chemically inert polymeric material

Cell Volume:

<1 μ L

Heat Exchanger:

Low dispersion

Maximum Cell Operating Pressure:

2 MPa (300 psi)

Suppressor Control:

AES, 0–150 mA in 1 mA increments

SRS, 0–500 mA in 1 mA increments

Column Heater

Operating Temperature Range:

Ambient +5 °C to 60 °C

Temperature Stability:

To ± 1 °C

Temperature Accuracy:

± 2 °C at 40 °C using external validation thermometer

Eluent Generator

Minimum and Maximum

Concentration:

0.1–100 mM

Flow rate:

0.1–3.0 mL/min

Maximum Operating Pressure:

21 MPa (3000 psi)

Maximum Solvent Concentration:

Anions: 25% methanol; Cations: no solvents

Physical Specifications

Power Requirements:

100–240 VAC, 50/60 Hz (power supply is autosensing, no voltage adjustment required)

Voltage Requirements:

90–265 VAC, 47–63 Hz

Operating Temperature:

4–40 °C (40–104 °F); cold-room-(4° C) compatible as long as system power remains on

Operating Humidity Range:

5–95% relative, noncondensing

Control Modes:

Full control through front panel and Chromeleon software; alternative control through TTL or relay closures; two relay out, two TTL out, four programmable inputs

USB Communication Protocol:

One input, built-in two-part USB hub

Leak Detection:

Built-in, optical sensor

Dimensions (h \times w \times d):

22.1 in \times 8.8 in \times 21 in

56.1 cm \times 22.4 cm \times 53.3 cm

Weight:

54 lbs

24.5 kg

Ordering Information

To order in the U.S., call (800) 346-6390 or contact the Dionex Regional Office nearest you. Outside the U.S., order through your local Dionex office or distributor. Refer to the following part numbers.

PART NUMBERS	
ICS-2000 Ion Chromatography System with Software and PC	
An ICS-2000/Chromeleon 6.5/Windows Workstation bundled package includes: an ICS-2000 with isocratic dual-piston pump, eluent generator to run Full EG, injection valve, column heater, heated conductivity cell, LCD touch-pad front panel, USB cable, Chromeleon version 6.50 PCS-2, Computer (with Windows XP or Windows 2000), and USB dongle. Comes with two Class 1 Timebases controlling one Dionex IC system. Consumables must be ordered separately.	
ICS-2000 Ion Chromatography System without	061087
Degas, Full EG, Chromeleon 6.5, Windows XP Workstation	
ICS-2000 Ion Chromatography System with	061086
Degas, Full EG, Chromeleon 6.5, Windows XP Workstation	
ICS-2000 Ion Chromatography System without	061089
Degas, Full EG, Chromeleon 6.5, Windows 2000 Workstation	
ICS-2000 Ion Chromatography System with	061088
Degas, Full EG, Chromeleon 6.5, Windows 2000 Workstation	



Reagent-Free is a trademark, and Atlas, AutoSuppression Chromeleon, and SRS are registered trademarks of Dionex Corporation. Windows is a registered trademark of Microsoft Corporation.



Dionex Corporation
 1228 Titan Way
 P.O. Box 3603
 Sunnyvale, CA
 94088-3603
 (408) 737-0700

Dionex Corporation
 Salt Lake City Technical Center
 1515 West 2200 South, Suite A
 Salt Lake City, UT
 84119-1484
 (801) 972-9292

Dionex U.S. Regional Offices
 Sunnyvale, CA (408) 737-8522
 Westmont, IL (630) 789-3660
 Houston, TX (281) 847-5652
 Atlanta, GA (770) 432-8100
 Marlton, NJ (856) 596-0609

Dionex International Subsidiaries
 Austria (01) 616 51 25 Belgium (32) 3-353 42 94 Canada (905) 844-9650 China (852) 2428 3282 Denmark (45) 36 36 90 90
 France 01 39 30 01 10 Germany 06126-991-0 Italy (06) 66 51 50 52 Japan (06) 6885-1213 The Netherlands (0161) 43 43 03
 Switzerland (062) 205 99 66 United Kingdom (01276) 691722
 * Designed, developed, and manufactured under an NSAI registered ISO 9001 Quality System.

